

## DIRECTIONS FOR USE:

One pack is designed for use by one patient and will enable them to cleanse their entire body. Nursing staff should assist if the patient is unable to use the wipes themselves. The wash cloths can be warmed or used straight from the pack. After the cleanse, allow your skin to air dry, no need to rinse off.

### Wash Mitts:



- 1.** For the 8 cloth packets, conduct full body cleansing as indicated in the diagram above.



- 2.** Clean the face, neck and hands.



- 3.** Use one wipe for the chest area. Use another wipe to clean the back.



- 4.** Use one wipe to clean the arm, hand and the armpit area. Repeat same for the other arm using a new wipe.



- 5.** Use one wipe to clean each leg.



- 6.** Use one wipe to clean groin, perineum and buttocks.

### Shampoo caps:



- 1.** Put on gloves prior using the product. There is one single use cap in a pack. Open the pack and place it over the patient's head.



- 2.** Massage it well and leave it for a few minutes for the product to cleanse and condition the hair.



- 3.** Remove it. No need to rinse. Discard used cap and gloves to bin.

## INGREDIENTS

**Wash Mitts:** Aqua, Dimethicone (3%), Linoleamidopropyl Pg-Dimonium Cholride Phosphate (Phospholipid), Lauryl Glucoside, Panthenol (Provitamin B5), Allantoin, Aloe Barbadensis Extract, Chamomilla Recutia Extract, Tocopheryl Acetate (Vitamin E), Chlorhexidine Digluconate, Glycerine, Polysorbate 20, Tetrasodium Edta, Citric Acid, Sodium Benzoate, Potassium Sorbate.

**Shampoo Caps:** Aqua, PEG-40 Glyceryl Cocoate, Phenoxyethanol, Sodium, Coceth Sulfate, Benzoate, Amodimethicone, Tetrasodium Glutamate, Diacetate, Parfum, Trideceth-12, Cetrimonium Chloride, Sodium Citrate