




MedVANCE[®]

Collagen Dressings



Improving patient care and enhancing the quality of people's lives



MedVANCE[®] Collagen Dressings



Product Description:

MedVance[®] Collagen Dressing Pads

MedVance Collagen Dressing Pads transform into soothing gel upon contact with exudate, absorbing wound drainage while maintaining an optimal moisture balance to aid in the formation of granulation tissue and epithelialization.

MedVance[®] Collagen Powder

MedVance Collagen Powder has a high surface area, allowing it to deliver the maximum amount of collagen to the wound surface. 100% natural collagen is easy to apply, absorbs light to moderate levels of wound exudate, and provides a moist wound healing environment.

Indications for Use:

MedVance[®] Collagen Dressings are intended for the management of low to moderately exuding acute or chronic, partial, or full thickness wounds including:

- Venous leg ulcers
- Diabetic Foot ulcers
- Pressure Injuries
- Second degree burns
- Superficial abrasions
- Donor sites
- Dehisced surgical wounds
- Traumatic wound healing by secondary intention

MedVance[®] Collagen Dressings

MedVance Collagen Dressings are 100% non-hydrolyzed type 1 bovine native collagen.

Proprietary technology process protects and retains native triple helical protein structure, allowing superior stability of the molecule and scaffolding through all four phases of wound healing.

MedVance Collagen Dressings are available in both powder and pad formats to fit a wide variety of wound types.

Features and Benefits:

- Intact (native) collagen provides highest bioactivity
- Compatible with topicals and secondary dressings
- Forms conformable gel on contact with fluid
- Supports moist wound healing environment
- Easy application and removal
- Shelf life of five years

Product Availability

MedVance[®] Collagen Dressing Pads

Product Code	Size
V7000202	2" x 2"
V7000404	4" x 4"

MedVance[®] Collagen Powder

V7000001	1g
----------	----

MedWay Inc.

1650 Horizon Parkway, Suite 450
 Buford, GA 30518
 (888) 5MEDWAY (633929)
 sales@medwayus.com
 www.medwayus.com

mv.br.7.1-v1

