




CURVANCE™

Foam Ag Dressing



Improving patient care and enhancing the quality of people's lives



CURAVANCE™ Foam Ag Dressing

Features and Benefits:

- Absorbent foam absorbs excess wound fluid to maintain the optimal moist wound healing environment
- Non-adherent foam reduces trauma during dressing change for enhanced patient experience
- Broad-spectrum antibacterial activity protects against the most common clinic pathogens associated with wound infection
- High MVTR semi-permeable PU film provides low coefficient of friction and extends wear time
- Transparent film serves as a physical barrier to prevent fluid and microbial contamination



Product Description:

CuraVance™ Foam Ag Dressings are multilayer, non-adherent, antibacterial foam dressings with silver, comprised of moisture vapor permeable waterproof backing, absorbent foam and a soft wound contact surface.

The semipermeable polyurethane film outer layer acts as a liquid and microbial barrier to prevent contamination while delivering balanced moisture vapor and oxygen transfer and a moist wound healing environment.

The non-adherent foam absorbs and retains exudate and ensures atraumatic dressing removal to deliver optimal patient outcomes.

Laboratory testing has demonstrated that CuraVance™ Foam Ag Dressings provide broad-spectrum antibacterial activity against the most common pathogens associated with wound infections for up to 7 days.


CuraVance™ Foam Ag Dressings are available with and without adhesive border.

Indications for use:

CuraVance™ Foam Ag Dressings are indicated for the management of moderate to heavily exuding, acute and chronic wounds, including:

- Diabetic foot ulcers
- Venous leg ulcers
- Pressure ulcers
- Graft and donor sites
- Traumatic wounds
- Second-degree burns
- Post-operative wounds and skin abrasions

Foam Ag - Sustained Antibacterial Activity



Product Availability

CuraVance™ Non-Bordered Foam Ag Dressing

Product Code	Size	Packaging	HCPCS
H1000404	4"x4"	5 pcs/box 80 box/case	A6209

CuraVance™ Bordered Foam Ag Dressing

H1010404	4"x4"	5 pcs/box 80 box/case	A6212
----------	-------	--------------------------	-------

MedWay Inc.

1650 Horizon Parkway, Suite 450
 Buford, GA 30518
 (888) 5MEDWAY (633929)
 sales@medwayus.com
 www.medwayus.com

cv.br.1.1-v1

