



EXTRA THICK WASH GLOVES



WEARABLE GLOVE - Just put your hand inside the glove and clean your body easily.
NO WATER, NO RINSE - Wash Gloves are designed for the waterless cleaning of the body and are very easy and fast to use without rinsing. They reduce contamination and are also heatable in the microwave if needed.



SPECIFICATIONS	DESCRIPTION
NONWOVEN	85 GSM WHITE HONEYCOMB SPUNLACE
SEWING:	Ultrasonic
PACKING:	FlowPack, 60 PE + 12 PET
FINISHED DIMENSIONS:	14,5cm x 22cm (Rectangular Shape)

PACKING DETAILS	PACK/BOX/PALLET AMOUNT	WEIGHT	DIMENSIONS
PACK	8 WASHGLOVES IN A PACK	0,24 KGS	
BOX	55 PACKS IN A BOX	14 KGS	39CM x 39 CM x 30CM
PALLET	36 BOXES IN A PALLET	520 KGS	80CM x 120CM x 195CM
DELIVERY TIME	20 FT CONTAINER / 14 DAYS		



SHAMPOO & CONDITIONER CAP

elivo Shampoo & Conditioner Cap is intended for cleaning the hair, head area, and scalp of bedridden and intensive care patients. With a maximum use of 3 minutes on the head, it ensures that maintenance and cleaning are completed in a short time. **To use**, first remove the cap from the pack, then place the cap on the head, making sure that all the hair is inside the cap, and begin massaging the hair from top to bottom for 3 minutes before removing the cap.

No rinse is required after using Shampoo Cap. It can be heated through microwave for 30 seconds at 700 W. The product contains many beneficial ingredients such as Vitamin E, Provitamin B5, Dimethicone, Allantoin, Glycerin, Chamomile Extract, Aloe Vera.



SPECIFICATIONS	DESCRIPTION
NONWOVEN	120 GSM WHITE SPUNLACE
OUTER NONWOVEN	PE FILM 35GSM WHITE-WATERPROOF
SEWING	Over-locked: 5 stitches per cm with T120 Natural Polyester
ELASTIC	PULAST Transparent Elastic, latex free, 300 mic. X 6mm
FINISHED DIMENSIONS	Component size: 50-53 cm diameter cut size
	Finished article: Aperture: 13-16cm diameter Diameter: 26-28cm diameter

PACKING DETAILS	PACK/BOX/PALLET AMOUNT	WEIGHT	DIMENSIONS
PACK	1 SHAMPOO CAP IN A PACK	0,160 KGS	50CM DIAMETER
BOX	75 PACKS IN A BOX	13 KGS	39CM x 39 CM x 30CM
PALLET	36 BOXES IN A PALLET	483 KGS	80CM x 120CM x 195CM

MedWay Inc.

1650 Horizon Parkway, Suite 450
 Buford, GA 30518
 (888) 5MEDWAY (633929)
 sales@medwayus.com
 www.medwayus.com