



**MedVANCE<sup>®</sup>**

*Super Absorbent Dressing*



Improving patient care and enhancing the quality of people's lives



# MedVANCE<sup>®</sup> Super Absorbent Dressing



## Product Description:

### MedVance<sup>®</sup> Super Absorbent Dressings

MedVance<sup>®</sup> Super Absorbent Wound Dressings are structured to absorb and retain significant amounts of wound exudate.

The hydrophilic wound contact layer does not adhere to the wound site while transmitting fluid into the dressing.

A distribution layer spreads exudate to maximize the use of the superabsorbent gelling fiber pad, enhancing fluid absorption and retention.

The semi-occlusive backing layer prevents fluid and microbial strike-through while maintaining an appropriate moist wound healing environment.

MedVance<sup>®</sup> Super Absorbent Wound Dressings are available with and without silicone adhesive border.

The silicone border acts as a gentle wound contact layer, which sticks to the skin around the wound but does not stick to the wound bed, providing atraumatic and pain-free removal at dressing changes.

## Features and Benefits:

- Gelling Fiber Core manages high volumes of exudate
- Locks in exudate under compression
- Highly conformable for patient comfort
- Distribution layer maximizes absorption capacity
- Support moist wound healing environment
- Easy application and removal

## Indications for Use:

MedVance<sup>®</sup> Super Absorbent Dressings are intended for highly exuding acute or chronic, partial, or full thickness wounds including:

- Venous leg ulcers
- Diabetic foot ulcers
- Pressure injuries
- Second-degree burns
- Superficial abrasions
- Donor sites
- Dehisced surgical wounds
- Traumatic wound healing by secondary intention

## Product Availability

MedVance <sup>®</sup> Super Absorbent Dressing			
Product Code	Dressing Size	Pad Size	Packaging
V6000303	3" x 3"	2.75" x 2.75"	5 pcs/box 40 box/case
V6000404	4" x 4"	3.75" x 3.75"	5 pcs/box 40 box/case
V6000606	6" x 6"	5.75" x 5.75"	5 pcs/box 40 box/case
MedVance <sup>®</sup> Bordered Super Absorbent Dressing			
V6010606	6" x 6"	4" x 4"	5 pcs/box 40 box/case

## MedWay Inc.

1650 Horizon Parkway, Suite 450  
 Buford, GA 30518  
 (888) 5MEDWAY (633929)  
 sales@medwayus.com  
 www.medwayus.com

mv.br.6.1-v1

