



# MEDVANCE<sup>®</sup> *Alginate*



Improving patient care and enhancing the quality of people's lives



# MedVANCE<sup>®</sup> Alginate



## Product Description:

MedVance<sup>®</sup> Alginate Calcium Alginate dressing is rich in soft gelling mannuronic acid.

The nonwoven pad format manages exuding wounds, and the rope is appropriate for packing cavities and indurations.

MedVance<sup>®</sup> Alginate dressings form a soft, hydrophilic gel immediately upon contact with wound exudate, to absorb excess exudate and inhibitors to wound healing and provide a moist wound healing environment in intimate contact with the wound bed.

Nonwoven reinforcing technology enhances wet strength, facilitating atraumatic and complete removal.

## Indications for use:

MedVance<sup>®</sup> Alginate dressings are indicated for the management of moderate to heavily exuding chronic or acute wounds, including:

- Venous leg ulcers
- Diabetic foot ulcers
- Pressure injuries
- Donor sites
- Partial and full thickness burns and abrasions

## Features and Benefits:

- Rapid gel formation—quickly absorbs excess exudate and inhibitors to wound healing
- Provides natural moist wound healing environment
- Intimately conforms to the wound bed, providing complete coverage
- Manages exudate under compression
- Extra nonwoven reinforcement ensures intact and complete removal
- Can be cut to shape and layered for best fit and increased absorption

Product Availability			
Product Code	Size	Packaging	HCPCS
V4000202	2"x2"	5 pcs/box 80 box/carton	A6196
V4000408	4"x8"	5 pcs/box 80 box/carton	A6197
V4004242	4.25"x4.25"	5 pcs/box 80 box/carton	A6197
V4000012	12" Rope	5 pcs/box 80 box/carton	A6199

MedWay Inc.

1650 Horizon Parkway, Suite 450  
Buford, GA 30518  
(888) 5MEDWAY (633929)  
sales@medwayus.com  
www.medwayus.com

mv.br.4.1-v1

