



# MEDVANCE<sup>®</sup>

## *Hydrocolloid*



Improving patient care and enhancing the quality of people's lives





### Features and Benefits:

- Provides a moist wound healing environment
- Supports autolytic debridement
- Comfortable to wear and easy to remove
- Protects the wound bed from external contamination and moisture-associated skin damage
- Low coefficient of friction reduces shear and skin stripping

### Product Description:

The MedVance® line of Hydrocolloid dressings includes Thin, Bordered, Foam-top products, as well as specialty shapes for difficult to dress locations.

Hydrocolloid dressings create a conformable gel on contact with wound exudate, providing a moist and insulating environment that supports wound healing and autolytic debridement.

The outer layer provides a barrier to fluid and microbial contamination, while the hydrocolloid layer adheres to periwound skin but not to the wound bed.

### Indications for use:

MedVance® Hydrocolloid dressings are indicated for the management of acute and chronic wounds presenting with low to moderate levels of exudate, including:

- Venous stasis ulcers (VLU)
- Pressure injuries
- Diabetic foot ulcers (DFU)
- Donor and graft sites
- Post-operative wounds
- Superficial and partial thickness burns, abrasions

### Product Availability

#### Hydrocolloid Thin Dressing

Product Code	Size	Packaging
V2400404	4"x4"	5 pcs/box 80 box/carton
V2400606	6"x6"	5 pcs/box 80 box/carton

#### Hydrocolloid Bordered Dressing

V2010202	2"x4"	5 pcs/box 80 box/carton
V2010404	4"x4"	5 pcs/box 80 box/carton
V2010606	6"x6"	5 pcs/box 80 box/carton
V2110607 Sacral	5"x7"	5 pcs/box 80 box/carton
V2213505 Heel	3.5"x5"	5 pcs/box 80 box/carton

#### Hydrocolloid Bordered Dressing with foam cushion

V2500404	4"x4"	5 pcs/box 80 box/carton
----------	-------	----------------------------

### MedWay Inc.

1650 Horizon Parkway, Suite 450  
 Buford, GA 30518  
 (888) 5MEDWAY (633929)  
 sales@medwayus.com  
 www.medwayus.com

mv.br.2.1-v1