



MEDVANCE[®]

Silicone Tape



Improving patient care and enhancing the quality of people's lives



MedVANCE[®] Silicone Tape

Features and Benefits:

- Highly conformable and extensible
- Comfortable to wear
- Unlike nonwoven tapes, extends to remain in place when stretched
- Extra-thick, medical-grade adhesive provides superior, reliable, durable adhesion
- Stays in place when others fail over joints, including heels, knees and elbows
- Reliable secondary and “window frame” seal-and-secure
- Superior holding power to secure tubing, cannulas and devices
- Transparent for sub-tape visualization of skin condition
- Low coefficient of friction to reduce friction, shear and skin tears
- Perforated for practical tear-to-length convenience
- Easy, atraumatic, pain-free removal without MARSI (Medical Adhesive Related Skin Injury)
- Silicone products may mitigate the appearance of scars during the remodeling phase
- Repositionable for optimal positioning and maximum utility under tubing and canulae



Product Description:

MedVance[®] Silicone Tape is designed to deliver enhanced conformability, extensibility, transparency and reliable adhesion through the use of extra-thick medical-grade silicone adhesive under a thin polyurethane backing.

MedVance Silicone Tape is waterproof, provides a low coefficient of friction, and tears easily for quick and practical application.

Premium hypoallergenic silicone adhesive remains securely in place while allowing atraumatic, pain-free removal with no adhesive residue.

Indications for use:

MedVance Silicone Tape is intended for use as a secondary dressing to secure and seal primary dressings, tubing, cannula and other medical devices, and to protect skin from such devices, as well as from friction, shear and skin tears.

Silicone has been employed to protect and reduce scars and to reduce MARSI (Medical Adhesive Related Skin Injury), pain and trauma on tape removal. Silicone adhesives allow for removal and repositioning of the tape.

Product Availability

MedVance[®] Silicone - Silicone Tape

Product Code	Size	Packaging	HCPCS
V5000115-1	1"x1.5 yards	1 roll/box 80 boxes/cs	A4452
V5000115-3	1"x1.5 yards	3 roll/box 80 boxes/cs	A4452
V5000150-1	1"x5 yards	1 roll/box 80 boxes/cs	A4452
V5000150-3	1"x5 yards	3 roll/box 80 boxes/cs	A4452
V5000215-1	2"x1.5 yards	1 roll/box 80 boxes/cs	A4452
V5000215-3	2"x1.5 yards	3 roll/box 80 boxes/cs	A4452
V5000250-1	2"x5 yards	1 roll/box 80 boxes/cs	A4452
V5000250-3	2"x5 yards	3 roll/box 80 boxes/cs	A4452

MedWay Inc.

1650 Horizon Parkway, Suite 450
 Buford, GA 30518
 (888) 5MEDWAY (633929)
 sales@medwayus.com
 www.medwayus.com

mv.ifu.5.1-v1

