



PRODUCT GUIDE



OUR BRANDS

MedVANCE[®] CURVANCE[™]





Introduction

It's our mission to provide innovative, quality wound care products with superior value to both patients and health care providers, improving patient care and enhancing the quality of people's lives.

We reduce your supply spending by delivering products that are quality-driven, so you can concentrate on what matters most: successful patient outcomes.

Partnering with Medway Group-USA Division:

- No Minimum Order Quantity
- Free Shipping
- Inventory Replacement Programs
- Drop Shipping Services - Direct to Patient (min. pkg fee)
- Increase Profit Margins with a Direct Account

Partnering with Medway Group-International Division:

For International Inquiry and Details, please contact:
sales@medwayus.com

Table of Contents

MedVance Advanced Wound Care

Alginates	5
Collagens.....	6
Foams	7
Hydrocolloids	8-9
Silicone Foams	10-11
Silicone Tape	12-13
Super Absorbents	14-15
Masks	16
Acne Patches.....	17

CuraVance Advanced Wound Care

Silver Alginate	19
Silver Foams.....	20-21
Silver Silicone Foams	22-23
SiloVance Nano Silver Gel & Wound Wash.....	24

Dymacare Hygiene Wipes

Bed Bath Mitts	26
Shower Caps	26

MedAir Mattress Systems

MedAir Elite	28
MedAir Comfort.....	29
MedAir Eco	29

Vakito NPWT Kits & Devices

NPWT Kits & Supplies	31
----------------------------	----

Contact Us.....	32
-----------------	----

Medvance



Product Description:

MedVance Calcium Alginate dressing is rich in soft gelling manuronic acid. The nonwoven pad format manages exuding wounds, and the rope is appropriate for packing cavities and indurations.

MedVance Alginate dressings form a soft, hydrophilic gel immediately upon contact with wound exudate, to absorb excess exudate and inhibitors to wound healing and provide a moist wound healing environment in intimate contact with the wound bed.

Nonwoven reinforcing technology enhances wet strength, facilitating atraumatic and complete removal.

Indications for Use:

MedVance Alginate dressings are indicated for the management of moderate to heavily exuding chronic or acute wounds, including: venous leg ulcers, diabetic foot ulcers, pressure injuries, donor sites, partial and full thickness burns and abrasions.

Product Availability:

Product Code	Size	Packaging	HCPCs
V4000202	2"x2"	5 pcs/box	A6196
V4000408	4"x8"	5 pcs/box	A6197
V4004242	4.25"x4.25"	5 pcs/box	A6197
V4000012	12" Rope	5 pcs/box	A6199



Product Description:

MedVance Collagen dressings are 100% non-hydrolyzed type 1 bovine native collagen. Proprietary technology protects and retains native triple helical protein structure, allowing for superior stability of the molecule and scaffolding through all four phases of wound healing.

MedVance Collagen dressings are available in both powder and pad formats to fit a wide variety of wound types.

Indications for Use:

MedVance Collagen dressings are intended for the management of low to moderately exuding acute or chronic, partial or full thickness wounds including: venous leg ulcers, diabetic foot ulcers, pressure injuries, second degree burns, superficial abrasions, donor sites, dehisced surgical wounds, traumatic wound healing by secondary intention.

Product Availability:

Product Code	Size	Packaging	HCPCs
V7000202	2"x2"	5 pcs/box	A6021
V7000404	4"x4"	5 pcs/box	A6021
V7000001	1 gram	5 pcs/box	A6010



Product Description:

MedVance Foam is a non-adherent polyurethane foam dressing that does not stick to the wound bed. The porous structure facilitates faster and stronger absorbent capacity and provides a moist wound healing environment. Breathable polyurethane film allows vapor transmission and extends wear time and the soft and comfortable foam may help redirect pressure.

MedVance Foam is available with and without hydrocolloid adhesive border.

Indications for Use:

MedVance Foam is intended for the management of moderate to heavily exuding acute and chronic wounds, including partial and full thickness wounds such as: diabetic foot ulcers, venous leg ulcers, pressure ulcers, donor sites, skin abrasions, superficial and partial thickness burns.

Product Availability:

MedVance Foam - Non Adherent Foam Dressing

Product Code	Size	Packaging	HCPCs
V1000202	2"x2"	5 pcs/box	A6209
V1000404	4"x4"	5 pcs/box	A6209
V1000606	6"x6"	5 pcs/box	A6210

MedVance Foam - Adhesive Bordered Foam Dressing

V1010404	4"x4"	5 pcs/box	A6212
V1010606	6"x6"	5 pcs/box	A6213



Product Description:

The MedVance line of Hydrocolloid dressings includes Thin, Bordered, Foam-Top, as well as specialty shapes for difficult to dress locations.

Hydrocolloid dressings create a conformable gel on contact with wound exudate, providing a moist and insulating environment that supports wound healing and autolytic debridement.

The outer layer provides a barrier to fluid and microbial contamination, while the hydrocolloid layer adheres to periwound skin but not the wound bed.

Indications for Use:

MedVance Hydrocolloid dressings are indicated for the management of acute and chronic wounds presenting with low to moderate levels of exudate, including: venous stasis ulcers, pressure injuries, diabetic foot ulcers, donor and graft sites, post-operative wounds, superficial and partial thickness burns, abrasions.

Product Availability:

Hydrocolloid Thin Dressing

Product Code	Size	Packaging	HCPCs
V2400404	4"x4"	5 pcs/box	A6234
V2400606	6"x6"	5 pcs/box	A6235

Hydrocolloid Bordered Dressing

V2010202	2"x2"	5 pcs/box	A6237
V2010404	4"x4"	5 pcs/box	A6237
V2010606	6"x6"	5 pcs/box	A6238
V2110607 Sacral	5"x7"	5 pcs/box	A6238
V2213505 Heel	3.5"x5"	5 pcs/box	A6237

Hydrocolloid Bordered Dressing with Foam Cushion

V2500404	4"x4"	5 pcs/box	A6212
----------	-------	-----------	-------



Product Description:

MedVance Silicone Foam dressing is an absorbent self-adhering dressing consisting of a polyurethane (PU) foam pad, a soft, hypoallergenic silicone adhesive wound contact layer, and a high moisture vapor transfer transparent film top layer.

The multilayer structure can absorb large amounts of wound exudate, and the soft silicone adhesive acts as a gentle wound contact layer that seals the wound edges but does not stick to the wound bed, minimizing pain and medical adhesive related skin damage (MARSI) when changing dressings.

Breathable PU film allows the passage of oxygen and water vapor, extending wear time and ensuring a moist wound healing environment.

MedVance Silicone Foam dressings are available with and without an adhesive border, and as a conformable exudate transfer dressing.

Indications for Use:

MedVance Silicone Foam Dressing is indicated for the management of moderate to highly exuding full or partial thickness wounds such as: pressure injuries, venous stasis ulcers, diabetic foot ulcers, graft or donor sites, abrasions,, and other wounds healing by secondary intention.

Product Availability:

MedVance Silicone Foam Dressing

Product Code	Size	Packaging	HCPCs
V3000202	2"x2"	5 pcs/box	A6209
V3000303	3"x3"	5 pcs/box	A6209
V3000404	4"x4"	5 pcs/box	A6209
V3000606	6"x6"	5 pcs/box	A6210
V3000405	4"x5"	5 pcs/box	A6210
V3000408	4"x8"	5 pcs/box	A6210
V3000608	6"x8"	5 pcs/box	A6210
V3000808	8"x8"	5 pcs/box	A6211

MedVance Silicone Foam Dressing with Border

V3010202	2"x2"	5 pcs/box	A6212
V3010205	2"x5"	5 pcs/box	A6212
V3010303	3"x3"	5 pcs/box	A6212
V3010404	4"x4"	5 pcs/box	A6212
V3010408	4"x8"	5 pcs/box	A6213
V3010408W Wings	4"x8"	5 pcs/box	A6213
V3010412	4"x12"	5 pcs/box	A6213
V3010412W Wings	4"x12"	5 pcs/box	A6213
V3010606	6"x6"	5 pcs/box	A6212
V3010608 Oval	6"x8"	5 pcs/box	A6213
V3010707	7"x7"	5 pcs/box	A6213
V3210508 Heel	5"x8"	5 pcs/box	A6212
V3210709 Heel	7"x9"	5 pcs/box	A6213
V3110707 Sacral	7"x7"	5 pcs/box	A6213
V3110909 Sacral	9"x9"	5 pcs/box	A6213

MedVance Silicone Foam Transfer Dressing

V3020608	6"x8"	5 pcs/box	A6210
V3020820	8"x20"	5 pcs/box	A6211

MEDVANCE[®] *Silicone Tape*



Product Description:

MedVance Silicone Tape is designed to deliver enhanced conformability, extensibility, transparency and reliable adhesion through the use of extra-thick medical-grade silicone adhesive under a thin polyurethane backing.

MedVance Silicone Tape is waterproof, provides a low coefficient of friction, and tears easily for quick and practical application.

Premium hypoallergenic silicone adhesive remains securely in place while allowing atraumatic, pain-free removal with no adhesive residue.

Indications for Use:

MedVance Silicone Tape is intended for use as a secondary dressing to secure and seal primary dressings, tubing, cannula and other medical devices, and to protect skin from such devices, as well as from friction, shear and skin tears.

Silicone has been employed to protect and reduce scars and to reduce MARSI (Medical Adhesive Related Skin Injury), pain and trauma on tape removal. Silicone adhesives allow for removal and repositioning of the tape.

Product Availability:

MedVance Silicone Tape

Product Code	Size	Packaging	HCPCs
V5000115-1	1"x1.5 yds	1 roll/box	A4452
V5000115-3	1"x1.5 yds	3 rolls/box	A4452
V5000150-1	1"x5 yds	1 roll/box	A4452
V5000150-3	1"x5 yds	3 rolls/box	A4452
V5000215-1	2"x1.5 yds	1 roll/box	A4452
V5000215-3	2"x1.5 yds	3 rolls/box	A4452
V5000250-1	2"x5 yds	1 roll/box	A4452
V5000250-3	2"x5 yds	3 rolls/box	A4452

MEDVANCE[®] *Super Absorbents*



Product Description:

MedVance Super Absorbent Wound Dressings are structured to absorb and retain significant amounts of wound exudate.

The hydrophilic wound contact layer does not adhere to the wound site while transmitting fluid into the dressing. A distribution layer spreads exudate to maximize the use of the superabsorbent gelling fiber pad, enhancing fluid absorption and retention. The semi-occlusive backing layer prevents fluid and microbial strike-through while maintaining an appropriate moist wound healing environment.

MedVance Super Absorbent Wound Dressings are available with and without silicone adhesive borders.

The silicone border acts as a gentle wound contact layer, which sticks to the skin around the wound but does not stick to the wound bed, providing atraumatic and pain-free removal at dressing changes.

Indications for Use:

MedVance Super Absorbent Dressings are intended for highly exuding acute or chronic, partial, or full thickness wounds including: venous leg ulcers, diabetic foot ulcers, pressure injuries, second-degree burns, superficial abrasions, donor sites, dehisced surgical wounds, traumatic wound healing by secondary intention.

Product Availability:

MedVance Super Absorbent Dressing

Product Code	Overall Size	Pad Size	Packaging	HCPCs
V6000303	3"x3"	2.75"x2.75"	5 pcs/box	A6196
V6000404	4"x4"	3.75"x3.75"	5 pcs/box	A6196
V6000606	6"x6"	5.75"x5.75"	5 pcs/box	A6197

MedVance Bordered Super Absorbent Dressing

V6010606	6"x6"	4"x4"	5 pcs/box	A6196
----------	-------	-------	-----------	-------



Product Description:

Disposable 3-layers face mask, Blue/White in color. The inner and outer layers are made of 100% spun-bond polypropylene. The middle layer is made of 100% melt blown polypropylene filter media. BFE (Bacterial Filtration Efficiency) = 95% which makes it easy to breath. Inner/outer Differential Pressure = 49 Pa/cm². Ear loops make it easy to wear. The nose piece contains a malleable aluminum wire to secure the mask while wearing. One size fits all.

Product Availability:

FDA510 Medical Grade Approved Masks

Product Code	Size	Packaging
MASK	7.5"x4"x3.5"	50 pcs/box

MedVANCE[®]

Acne Patches



Product Description:

Stubborn acne or a surprise breakout – no problem! MedVance Acne Cover Patches actively work to heal blemishes and shrink acne spots away. The hydrocolloid absorbs the pimple fluid within 6-8 hours of application helping to heal the blemish quickly and discretely. Eliminates picking and popping and leaves a clearer looking skin. The thin profile of the patch blends well with all skin tones. Clinically tested, drug-free and never tested on animals.

**Product Availability:
MedVance Acne Cover Patches**

Product Code	Size	Packaging
V8000814-4	multiple	144 patches/pack

CURAVANCE



Product Description:

CuraVance Silver Alginate antibacterial calcium alginate dressing is rich in highly absorbent soft gelling mannuronic acid. The non-woven pad format manages exuding wounds, and the rope is appropriate for filling cavities and indurations. CuraVance Silver Alginate immediately forms a soft, hydrophilic gel on contact with wound exudate absorbing excess fluid and exudate-associated inhibitors of wound healing, creating a moist wound healing environment in intimate contact with the wound bed. Non-woven reinforcing technology enhances wet strength, facilitating atraumatic and complete removal. CuraVance Silver Alginate provides antibacterial protection against both gram-negative and gram-positive bacteria, commonly associated with wound infections.

Indications for Use:

CuraVance Silver Alginate is indicated for the management of moderate to heavily exuding chronic or acute wounds, including venous leg ulcers, diabetic foot ulcers, pressure injuries, donor sites, partial and full thickness burns and abrasions.

Product Availability:

Product Code	Size	Packaging	HCPCs
H4000202	2"x2"	5 pcs/box	A6196
H4004242	4.25"x4.25"	5 pcs/box	A6197
H4000012	12" Rope	5 pcs/box	A6199

CURAVANCETM Silver Foams



Product Description:

CuraVance Foam Ag Dressings are multilayer, non-adherent anti-bacterial foam dressings with silver, comprised of moisture vapor permeable waterproof backing, absorbent foam and a soft wound contact surface. The semipermeable polyurethane film outer layer acts as a liquid and microbial barrier to prevent contamination while delivering balanced moisture vapor and oxygen transfer and a moist wound healing environment. The non-adherent foam absorbs and retains exudate and ensures atraumatic dressing removal to deliver optimal patient outcomes.

CuraVance Foam Ag Dressings provide broad-spectrum antibacterial activity against the most common pathogens associated with wound infections for up to 7 days.

CuraVance Foam Ag Dressings are available with and without adhesive border.

Indications for Use:

CuraVance Foam Ag Dressings are indicated for the management of moderate to heavily exuding, acute and chronic wounds, including: diabetic foot ulcers, venous leg ulcers, pressure ulcers, graft and donor sites, traumatic wounds, second-degree burns, post-operative wounds and skin abrasions.

Product Availability:

CuraVance Non-Bordered Foam Ag Dressing

Product Code	Size	Packaging	HCPCs
H1000404	4"x4"	5 pcs/box	A6209

CuraVance Bordered Foam Ag Dressing

H1010404	4"x4"	5 pcs/box	A6212
----------	-------	-----------	-------

CURAVANCETM Silver Silicone Foams



Product Description:

CuraVance Silicone Ag Foam Dressings are multilayer antibacterial foam dressings incorporating a moisture vapor permeable water-proof backing, absorbent Ag foam, a super-absorbent core and soft silicone skin, and a wound contact surface. Dressing wear time is enhanced through the use of a super-absorbent layer and a fluid transfer layer that ensures maximum use of absorbent capacity and retention of fluid under compression.

A semipermeable polyurethane film outer layer acts as a liquid and microbial barrier to prevent both strike-through and contamination while delivering balanced moist wound healing. Soft, silver-impregnated foam absorbs and retains exudate while providing broad-spectrum antibacterial activity for up to 7 days. The soft silicone adhesive wound contact layer and border maintains patient comfort, avoids MARSI, and provides pain-free atraumatic removal.

Indications for Use:

CuraVance Silicone Ag Foam Dressings are intended for the management of partial and full thickness, acute and chronic wounds with moderate to heavy exudate, including: diabetic foot ulcers, venous leg ulcers, pressure injuries, second-degree burns, skin graft and donor sites, surgical and traumatic wounds, post-operative wounds and skin abrasions.

Product Availability:
CuraVance Silicone Foam Ag Dressing

Product Code	Size	Packaging	HCPCs
H3010404	4"x4"	5 pcs/box	A6212
H3110707	7"x7"	5 pcs/box	A6213
H3110909	9"x9"	5 pcs/box	A6214

CURA**VANCE**TM Silver Gel & Wash



Product Description:

SiloVance Antibacterial Wound Dressing Gel is a water-based gel wound dressing for use in topical OTC management of: first and second degree burns, lacerations, abrasions, minor cuts, skin irritations, and by or on the order of a licensed healthcare practitioner for: wounds such as stasis ulcers, pressure ulcers, diabetic ulcers, surgical incision sites, device insertion site wounds, skin tears, graft sites and donor sites. SiloVance Antibacterial Wound Gel contains silver that in laboratory tests has been shown to inhibit the growth of microorganisms such as Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, antibiotic-resistant bacteria, MRSA and VRE as well as fungi such as Candida albicans.

Product Availability:

SiloVance Nano Silver Gel

Product Code	Size	Packaging	HCPCs
H5200001	1oz	1 pcs/box	A6260
H5200003	3oz	1 pcs/box	A6260

SiloVance Nano Silver Wound Wash

Product Code	Size	Packaging	HCPCs
H5300004	4oz	1 pcs/box	A6248

The logo features the word "eLivo" in a white, bold, sans-serif font. A small white teardrop icon is positioned above the final 'o'. A thin white curved line, resembling a smile, is located directly beneath the text. The entire logo is centered within a dark blue, irregular, blob-like shape. This central shape is layered over a larger, lighter blue, similarly irregular shape, creating a sense of depth. The background is a solid dark blue.

eLivo



Hygiene Wipes



Wash Mitts

Fragranced

Lightly scented skin washing mitt, enriched with Aloe Vera. Suitable for all skin types. Excellent for bed bath.

#E2112628 / 3-packs of 8 gloves ea. (24)



Shampoo Cap

Fragranced

Cleanses and cares for the hair and scalp. A moist mild shampoo and conditioner. Suitable for all hair types. No need to rinse, simply towel dry.

#E2111423 / 3-packs of 1 cap ea. (3)



MedAir®

The logo features a stylized graphic above the text. The graphic consists of a large, open circle at the top, with two wavy lines flowing from its base. A smaller, tighter swirl is positioned to the right of the main wavy lines, just before the registered trademark symbol. The text "MedAir" is written in a white, italicized, sans-serif font.



Mattress Systems



MedAir Elite

The **MedAir Elite** Dual Therapy Group II mattress system delivers immersive alternating pressure and low air loss micro-climate management in a versatile and cost-effective design. The fully customizable digital pump adjusts air flow and pressure redistribution to individual patient requirements, while diagnostic alarm settings inform caregivers when patient attention is required. Auto return function automatically shifts from Static to Alternating Pressure mode to ensure maximum pressure redistribution and relief. The 8" deep Cell-on- Cell configuration with Pillow Function provide immersive pressure relief and prevents bottoming out, while providing 24 hr. power outage protection.

Adjustable low air loss settings deliver low interface pressures and regulate temperature and moisture microclimate to reduce the most common factors leading to pressure injuries and the impediments to healing, optimizing patient comfort, wound healing and pressure injury prevention.

Product code: MA-M4P8-US

Specifications

Dimensions: 80"(L) x 36"(W) x 8"(H);

Pump Airflow: 8 LPM;

Cycle Time: 10, 15, 20, 25 Minutes

Weight Capacity: 450 lbs.



Mattress Systems



MedAir Comfort

MedAir Comfort mattress system is 4.5" deep with an 18 cell-on-cell configuration. The adjustable pump adjusts air flow and pressure redistribution.

#MA-CMFRT

Size: 78"Lx36"Wx4.5"D

Weight capacity: 360 lbs.



MedAir Eco

MedAir Eco mattress bubbles are on A-B lines which alternate inflating to provide superior therapeutic treatment. Made of durable heavy low toxic premium vinyl.

#MA-ECO

Size: 74"Lx36"Wx2.5"D

Weight capacity: 300 lbs.

VAKITO  **TM**



NPWT Dressing Kits

VAKITO NPWT Kits

Universal NPWT kits are available in Small, Medium and Large kit sizes.

Each kit contains individually wrapped and sterilized foam, drape and suction bell.

All Components are available for individual sale

VAKITO Kits are a cost-effective alternative



High quality products to help heal your patient's wounds!

NPWT Kits	
Item#	Description
MVF-KIT-S	Small NPWT Kit - 1 Suction Bell, 1 Drape, 1 Foam (20 kits per case)
MVF-KIT-M	Medium NPWT Kit - 1 Suction Bell, 1 Drape, 1 Foam (20 kits per case)
MVF-KIT-L	Large NPWT Kit - 1 Suction Bell, 1 Drape, 1 Foam (20 kits per case)
MVF-FOAM-S	Small Foam Single Pack Dressing (10x7.5x2.5cm)
MVF-FOAM-M	Medium Foam Single Pack Dressing (18x12.5x2.5cm)
MVF-FOAM-L	Large Foam Single Pack Dressing (26x15x2.5cm)
MVF-DRAPE	Drape Single Pack (30x30cm)
MVF-BELL	Suction Bell Single Pack (Round port with tube 1.5M)
MV-CON-M	Connection Tube Single Pack
MV-SIL-3015	Silicone Contact Sheet (30 x 15cm - 12"x6")



Contact Us

Mailing & Billing Address:

Medway Inc.
1650 Horizon Parkway, Suite 450
Buford, GA 30518

Customer Service: 888-5MEDWAY (633929)

USA Sales Contacts:

Andrea Nicolai, Sr. Business Development Manager
anicolai@medwayus.com
404-822-7860

Michele Endreszl
Independent Representative
Reserve Medical, llc
mendreszl@medwayus.com
sales@reservemed.com
470-894-9657 / 248-431-4054

